

## Apolipoprotein B Assay Kit (APOB)

**Method:** Immunoturbidimetric

### INTENDED USE

Cat . No.	Size	Instrument
GB170Z	R1: 4×45 ml R2: 4×15 ml	For Hitachi 717 & ShimadzuCL7200/8000
GS171Z	R1: 3×60 ml R2: 3×20 ml	For Hitachi917 & OlympusAU640/400/600
GX171Z	R1: 2×60 ml R2: 2×20 ml	For SYNCHRON CX4-5-7-9/LX20/DXC600-800

For the *in vitro* quantitative determination of apolipoprotein B (Apo B) in serum. For use as an aid in assessing the risk of cardiovascular disease.

### CLINICAL SIGNIFICANCE

Lipids are synthesised in the intestine or liver but must be transported to tissues and organs. However this is not possible without hydrophilic adaptation. Lipids are therefore transported by a series of micellar structures. These structures consist of an outer monolayer of protein (an apolipoprotein) and polar lipids (phospholipids and unesterified cholesterol) plus an inner core of neutral lipids (triglycerides and cholesterol esters). The apolipoproteins interact with a series of enzymes and tissue receptors and are therefore responsible for further metabolism and catabolism of the micelle.

The B apolipoproteins are the main form of protein found in low density lipoproteins (LDL). Two forms of apo B are found in humans. The most common form is apo B-100 (or large B) which constitutes the apo B found in lipoproteins synthesised in the liver. The other form is apo B-48 (or small B), thought to be synthesised in the intestinal wall. Apo B is the main cholesterol carrying protein in the blood and is the ligand concerned with the uptake of cholesterol into cells by the LDL-receptor pathway. Apo B shows atherogenic signs and is thus useful for the evaluation of coronary risk. Studies have shown that there is an inverse relationship between APO A-1 and coronary artery disease and a direct relationship with APO B such that patients with CAD have generally reduced levels of APO A-1 and increased levels of APO B.

### ASSAY PRINCIPLE

This method is based on the reaction of a sample containing human apo B and a specific antiserum to form an insoluble complex which can be measured turbidimetrically at 340 nm. By constructing a standard curve from the absorbances of standards the concentration of apo B can be determined.

### REAGENT COMPOSITION

Contents	Concentrations of reagents
<b>R1. Buffer</b>	
Polyethylene glycol	maximum 4%
Tris/HCl buffer	17 mmol/L
Sodium Chloride	125 mmol/L

<b>R2. Antibody Reagent</b>	
Anti-human-apo B	
<b>Calibrator</b>	lot specific

### SPECIMEN COLLECTION AND PREPARATION

Serum samples.


Apo B is stable for 7 days at 15 -25°C, 4 weeks at 2-8°C or 2 months at - 20°C.

### STABILITY AND PREPARATION OF REAGENTS

Reagents are stable until expiry date when stored unopened at 2-8°C.

### ASSAY PROCEDURE

Wavelength(main/sub) : 340 nm/700 nm



### CALIBRATION

Apolipoprotein Calibrator (value is lot specific) provided with the kit is recommended for calibration.

Recalibration is recommended for each series of samples run. The Apolipoprotein B concentration was determined using a WHO/IFCC reference standard.

### CALCULATIONS OF RESULTS

Plot calibrator concentrations against the corresponding-A values using graph paper. The concentration of Apo B in the sample is obtained by reading off -A value from the calibration curve. Exploration above or below the range of calibrators, must not be attempted.

### QUALITY CONTROL

Specific Protein Controls are recommended. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error, the following steps should be taken:

1. Check instrument settings and light source.
2. Check reaction temperature.
3. Check expiration date of kit and contents.

### NORMAL VALUE

Adults: 63 - 114 mg/dl

It is recommended that each laboratory should assign its own normal range as this is dependent upon geographical location.

### MAIN PERFORMANCE CHARACTERISTICS

#### LINEARITY

The method is linear up to 230 mg/dl. If the concentration is above 230mg/dl, please dilute it with 0.9% NaCl and reassay. Multiply the result by dilution factor.

Beijing Strong Biotechnologies, Inc.

Add: 5/F Kuang Yi Building, No. 15 Hua Yuan Dong Lu, Haidian District, Beijing 100191 P. R. China

Tel: +86 10 8201 2486 Fax: +86 10 8201 2812

Web: [www.bsbe.com.cn](http://www.bsbe.com.cn) Email: [jq.tech@bsbe.com.cn](mailto:jq.tech@bsbe.com.cn)

## PRECISION

The CV of the test should be less than 5%.

Intra assay precision		
N=20	level 1	level 2
Mean(mg/dl)	59.35	61.11
SD	0.84	1.00
CV(%)	1.42	1.64

Inter assay precision			
N=5	Batch 1	Batch 2	Batch 3
Mean(mg/dl)	57.45	59.01	59.20
$\bar{x}$	58.55		
$(X_{max}-X_{min})/\bar{x}$	$(59.2-57.45)/58.55 \times 100 = 2.99\%$		

## INTERFERENCE

The following analytes were tested up to the levels indicated and found not to interfere:

Hemoglobin	up to 500 mg/dl
Intralipid	up to 500 mg/dl
Bilirubin	up to 40 mg/dl

## SAFETY PRECAUTIONS AND WARNINGS

- For in vitro diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
- Reagent contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
- Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide from building up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
- Specimens should be treated as potentially infectious (HIV, Hepatitis B virus, Hepatitis C virus, etc.) and handled with appropriate caution.
- Reagents with different lot numbers should not be interchanged or mixed.

## REFERENCES

- Adolphson, J.L., Albers J.J., Journal of Lipid Research 1989; 30:599-806
- Bachorik, P.S., Kwiterovich, P.O. Clinica Chimica Acta 1988; 178:1-34
- Labour, C., Shepherd, J., Rosseneu, M. Clin. Chem. 1990; 36(4):591-597
- Wang, X.L., Dudman, N.P.B. Clin. Chem. 1989; 35(10):2082-2086
- Wang, X.L., Dudman, N.P.B., Wilken, D.E.L., Clin Chem 1989; 35(6):1000-10004
- Marcinova, S.M. et al (1992), WO/IFCC Meeting on Standardisation of Apolipoproteins, May 1992, Nice,

France

- Albers, J.J. et al (1992) Clin Chem. 38:658
- Provisional normal values recalculated on the basis of the CDC values C. Fruchart, J-C. (1986), Ann. Biol. Clin. 44:116
- Kukita H., Hiwada K., Kobubu T. Serum Apolipoprotein A-1, A-2 and B levels and their discriminative values in relatives of patients with coronary artery disease. Atheroscler 51:261, 1984.

## INDEX OF SYMBOLS



Manufacture



Catalogue Number



Lot number



Date of manufacture



Use by (Expiration date)




For In-Vitro Diagnostic use only



Stored at 2-8°C



Attention: See instruction for use



Authorized Representative in the European Company

Manufacture: Beijing Strong Biotechnology, Inc.

Address: No. 15, Yanqi North Second Street, Yanqi Economic Development Area, Huairou District, Beijing 101400, P. R. China

Tel: +86 10 61667168

EC REP : Lotus NL B.V.

Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

E-mail: peter@lotusnl.com

Tel: +31645171879 (English), +31626669008 (Dutch)

Beijing Strong Biotechnologies, Inc.

Add: 5/F Kuang Yi Building, No. 15 Hua Yuan Dong Lu, Haidian District, Beijing 100191 P. R. China

Tel: +86 10 8201 2486 Fax: +86 10 8201 2812

Web: [www.bsbe.com.cn](http://www.bsbe.com.cn) Email: [jq.tech@bsbe.com.cn](mailto:jq.tech@bsbe.com.cn)



Revised 2014 version 0801